

Flujómetros Integrados de Oxígeno y de Aire medicinal


Características

El Flujómetro Integrado Amvex combina un Flujómetro y una Toma de gas medicinal en un diseño único y compacto

- Perilla para control de flujo ajustable y fácil de usar
- Igual al usado tradicionalmente, solo que éste no permite ser removido
- El Etiquetado de código de color en la toma y en el flujómetro hace que sea rápido y fácil de distinguir entre Oxígeno y el Aire medicinal
- Disponible en diferentes tasas de flujo, tipo de gas y en tres diferentes estilos de conectores.



FI-T15UO-ON-E



DIAL STYLE INTEGRATED FLOWMETER CONFIGURATOR:

FI-TXXYZ-WN-S

Flowrate:

200 CCM:	2C
1 LPM:	01
3.5 LPM:	03
8 LPM:	08
15 LPM:	15

Note: Medical Air is only available in 0-15 LPM

Color Code:

USA:	U
ISO:	I

Gas Type:

Oxygen:	O
Medical Air:	A

Line Fitting:

DISS Male with demand check valve:	D
Ohmeda Coupler:	O
Chemetron Coupler:	C

Outlet Compatibility:

Diamond 1*:	A
Diamond 2&3 Beacon Medaes Diamond Care**:	B
Ohio Medical® Amico Beacon Series B Gentec Oxequip/Medstar Series #7 or 5/8-24*:	E
Puritan Bennett Series 73* Puritan Bennett Series "A"* Puritan Bennett Litharge* NCG/Chemetron Series 500, 560* Oxequip/Medstar Smooth Bore or 5/8-27*:	F
NCG/Chemetron Series 148, 400, 460*:	G
Tri-Tech/Ohmeda:	T